

Dairy Connection, Inc. 2312 Vondron Rd., Madison WI 53718 PH: (608) 242-9030 Fax: (608) 242-9036 Email: info@dairyconnection.com Web: dairyconnection.com Web: dairyconnection.com

PRODUCT DESCRIPTION SHEET NEC-C3 FRESH CHEESE STARTER CULTURE 1 DOSE / CH-3225-DC

Updated 05.12.23

DESCRIPTION

NEC-C3 is a blend of undefined, mesophilic multiple-species freeze-dried strains of DL-lactic acid bacterial cultures for making Fresh Cheese.

USAGE AND HANDLING

Products made from this using recommended dosage levels.

DIRECTIONS FOR USE

Shake well before use. Add directly to the milk avoiding foam and air introduction in the milk.

STORAGE AND SHIPPING

Use within 12 months from date of packaging only when stored at or below freezing. "Best used by" date is printed on label for reference.

Store in a dry atmosphere. Check before use that culture is in free-flowing powder form. If the product has formed a solid mass, it should be discarded.

Ambient temperatures during a 3-5 day transit time is acceptable.

COMPOSITION

Lactococcus lactis subsp. lactis Lactococcus lactis subsp. cremoris Lactococcus lactis subsp. lactis biovar. diacetylactis Leuconostoc mesenteroides subsp. cremoris Carrier: lactose

PROPERTIES

- Mesophilic homofermentative culture
- Mesophilic aroma culture developing diacetyl flavor and CO₂
- Direct vat inoculation
- Standardized activity

MICROBIOLOGICAL SPECIFICATIONS

Microbiological quality control - standard values.

Yeasts and Molds	<50 CFU/g
Coagulase – positive Staphylococci	<10 CFU/g
Salmonella sp.	neg. /25g
Listeria sp.	neg. /25g

PACKAGING

Product is packaged in layered "Biaxially oriented polypropylene" (BOPP) and white "Metallized" (MET-OPP) film.

PURITY AND LEGAL STATUS

Local regulations should always be consulted concerning the status of this product, as legislation regarding its intended use may vary.

SAFETY AND HANDLING

Safety Data Sheet (SDS) is available on request.

COUNTRY OF ORIGIN USA

GMO STATUS

To the best of our knowledge and according to documentation received from our suppliers, NEC-C3 does not contain any ingredients that contain genetically modified organisms.

The information contained in this publication is based on information provided by the manufacturer's research and development work and is to the best of our knowledge reliable, per the manufacturer of this product. Users should, however, conduct their own tests to determine the suitability of this product for their own specific purposes and the legal status for their intended use of the product. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for the infringement of any patents. Page 1 of 2



Dairy Connection, Inc. 2312 Vondron Rd., Madison WI 53718 PH: (608) 242-9030 Fax: (608) 242-9036 Email: info@dairyconnection.com Web: dairyconnection.com

PRODUCT DESCRIPTION SHEET NEC-C3 FRESH CHEESE STARTER CULTURE 1 DOSE / CH-3225-DC Updated 05.12.23

ALLERGENS

Table below indicates the presence (as added components) of the following allergens and products thereof:

YES	NO	ALLERGENS	DESCRIPTION OF COMPONENTS
	Х	Wheat	
	Х	Other cereals containing	
	Х	gluten Crustacean shellfish	
	Х	Eggs	
	Х	Fish	
	Х	Peanuts	
	Х	Soybeans	
х		Milk (including lactose)	
	Х	Nuts	
	Х	Celery	
	Х	Mustard	
	Х	Sesame seeds	
	Х	Sulphur dioxide and sulfites (>10 mg/kg)	
	Х	Lupin	
	Х	Mollusks	
Local regulation has always to be consulted as allergen labeling requirements may vary from country to country.			

The information contained in this publication is based on information provided by the manufacturer's research and development work and is to the best of our knowledge reliable, per the manufacturer of this product. Users should, however, conduct their own tests to determine the suitability of this product for their own specific purposes and the legal status for their intended use of the product. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for the infringement of any patents. Page 2 of 2